



Elim Christian School
VAGAL NERVE STIMULATOR CARE PLAN

Directions for magnet use

(Updated: 07/31/2015)

Name _____ Date _____

My magnet is kept: _____

My VNS should be swiped when:

Swipe the magnet across the generator located _____

_____, under the skin, slowly counting to _____
(one-one hundred, two-one hundred, three-one hundred etc.)

Do not hold the magnet over the generator for longer than 5 seconds.

May repeat _____ times.

Wait _____ seconds between swipes

Additional Instructions: _____

Keep "swiping magnet" away from items that are affected by magnetic fields like: desk top computers, computer disks, lab tops/iPads, credit cards, and batteries.

Keep any item that has a magnet at least 6 inches away from implanted device. Any magnet can activate and/or deactivate the device; i.e. ipads, speakers etc. Keep these items at desk or lap level.

Parent Signature

Printed Name